



# MASTER LEDtube EM/ Mains

## MASTER LEDtube Starter EMP GenII OP

The Philips MASTER LEDtube integrates a LED light source into a traditional fluorescent form factor. Its unique design creates a perfectly uniform visual appearance which cannot be distinguished from traditional fluorescent. For those that are looking for value for money within limited budget and re-lamping efforts for better light effect and lifetime.

### Warnings and Safety

• -

### Product data

General Information	
EU RoHS compliant	Yes
Product Data	
Full product code	871869681965400
Order product name	MASTER LEDtube Starter EMP GenII OP
EAN/UPC – product	8718696819654

Order code	81965400
SAP numerator – quantity per pack	1
Numerator – packs per outer box	40
SAP material	929001801432
SAP net weight (piece)	0.004 kg

# MASTER LEDtube EM/Mains

## Dimensional drawing



Product	D1	D2	C1
MASTER LEDtube Starter EMP GenII OP	3 mm	21.5 mm	34.5 mm

Dummy starter for cost down version

